


caBIG/ClinTrials Workspace/Lab Interface SIG

Lab Testing Ontology: Learning from HL7

Donald Connelly, MD, PhD, FCAP, FASCP, FACMI
University of Minnesota Cancer Center
September, 2004





Overview



- Purpose: what can we learn from HL7 about handling laboratory test results?
- Introductory lab testing terms
- HL7 v2.x Laboratory result messages
- What about lab results & LOINC?
- What will HL7 v3 do to what we do now?



Introductory lab testing terms

- Orderable Test (e.g., Blood Gas)
 - Appears on lab order forms
 - has a local order code (e.g., ABG for Arterial Blood Gas)
 - has a CPT (Current Procedural Terminology) code
- Non-orderable Test (e.g, pH of arterial blood)
 - A component test of an orderable test
 - PaO_2 , PaCO_2 , O_2 , SaO_2 , HCO_3 and pH are the components of the orderable "Blood Gas"
 - Sodium, Potassium, Bicarb, Chloride are the non-orderable components of the "Electrolytes" orderable test
 - usually doesn't have a local order code so can't be ordered by itself
 - reportable even though not orderable
 - reported with other components of an orderable test
 - Non-orderable components often dictated by lab analyzer technology and/or local practice
- Battery, Panel, Profile (e.g., Admission Battery)
 - An orderable test encompassing a group of orderable tests
 - An Admission Battery might be comprised of the four following orderable tests; CBC, Lytes, AST, and Urinalysis

University Hospitals**Laboratory Services**

11460*13108 Euclid Avenue • Cleveland, OH 44106-5199

(216)844-5227 • FAX (216)844-7560

RUBEN, BILLY

MRN: 3342233

BIRTHDATE SEX: 12/28/19 F

PATIENT #: 02351804

ORDERED: 07/16/01 10:56

ROOM# BED:

**This has been
HIPAA-fied**REQUESTED BY: OVERMOYER, BETH,
ATT. PHYSICIAN: OVERMOYER, BETH,

ADDT'L COPY TO:

FINAL

TEST NAME	RESULTS	REF VALUES	UNITS	TS
-----------	---------	------------	-------	----

HEMATOLOGY

COLLECTED 07/16/01 10:56

COMPLETE BLOOD COUNT


WBC	H	16.4	4.4-11.3	X10E9/L	UH
RBC		4.02	4.00-5.20	X10E12/L	UH
HGB	L	10.6	12.0-16.0	g/dL	UH
HCT	L	33.0	36.0-46.0	%	UH
MCV		82	80-100	fL	UH
MCH		26.4	26.0-34.0	pg	UH
MCHC		32.1	32.0-36.0	%	UH
PLT		228	150-450	X10E9/L	UH
RDW-CV	H	15.8	11.5-14.5	%	UH
MPV		10.6	7.2-11.1	fL	UH

DIFFERENTIAL

SEG NEUTROPHIL		21.0		%	UH
BAND NEUTROPHIL		28.0		%	UH
LYMPHOCYTE		11.0		%	UH
MONOCYTE		17.0		%	UH
EOSINOPHIL		0.0		%	UH
BASOPHIL		0.0		%	UH
METAMYELOCYTE		10.0		%	UH
MYELOCYTE		13.0		%	UH
SEG NEUTROPHIL		3.44	1.20-7.00	X10E9/L	UH
BAND NEUTROPHIL	H	4.59	0.00-0.70	X10E9/L	UH
LYMPHOCYTE		1.80	1.20-4.80	X10E9/L	UH
MONOCYTE	H	2.79	0.10-1.00	X10E9/L	UH
EOSINOPHIL		0.00	0.00-0.70	X10E9/L	UH
BASOPHIL		0.00	0.00-0.10	X10E9/L	UH
METAMYELOCYTE		1.64	Occ	X10E9/L	UH
MYELOCYTE		2.13	Occ	X10E9/L	UH

continued on next page

**Ordered
Test
Codes****CBC****DIFF**



What is HL7? (Health Level Seven)

- an application protocol for data exchange between computer systems in health care environments
- defines the messages as they are exchanged and the procedures used to exchange them
- independent of operating systems and programming languages



HL7 messages are made up of segments

- HL7's Unsolicited Observation Message (ORU) is sent out by the LIS when a lab result is completed
- The ORU Message must contain the following segments:
 - MSH Message Header (information about the message)
 - PID Patient Identification (information about the patient)
 - OBR Observation Request (information about the order)
 - OBX Observation Report (information about a result)
- Sometimes the ORU will also have PV (patient visit), ORC (order common) and NTE (note) segments but lets keep things simple

Message Subject Heading (MSH) Segment

- All messages begin with a MSH segment
 - Defines the intent, source, destination and syntax of a message
- Contains fields such as
 - Sending Application (MSH.3)
 - Sending Facility (MSH.4)
 - Receiving Facility (MSH.6)
 - Message Type (MSH.9)
- MSH example
- MSH|^~\&|PATHLAB|MERCY||GP009|200207251030||ORU^R01|CNTRL-3458|P|2.4

HL7 Tip

Fields are separated by vertical bars “|”

The Patient Identification Segment (PID)

- PID segment includes;
- includes fields such as:
 - Patient Identifier list (PID.3)
 - Patient Name (PID.5)
 - Date/Time of Birth (PID.7)
 - Administrative Sex (PID.8)
 - Patient Address (PID.11)

- PID example

PID|||666442222^^^MERCY^MR||MULLINS^LARRY||1956
0312|M|||6 ASHWOOD DRIVE^^RATHFARNHAM^DUBLIN

HL7 Tip

Fields may consist of components separated by ^

For example, patient name field consists of
|LastName ^ FName ^

MI

ORU Laboratory Report Message

- MSH| ...
 - PID| ...
 - OBR|1||1045814^PATHLAB|85022^HB|||200207250830||||
||||75555^COLLINS^PATRICIA^^^MB|||||||F
- Ordered Test Collection Date/Time
- Observation identifier
- Ordering Provider
- OBX|1|ST|718-7^HEMOGLOBIN:^LN||13.4|GM/DL|14-18|
N||S|F
- Abnormal Flag Result Status (F is final) Value Units Ref Range








Blood Gas example

(with a note segment, **NTE**)

- MSH ...
- PID ...
- OBR|1||1045825^LAB|L2560.000^BLOOD GAS, ARTERIAL (R)|||19940320105800|""|
1^EA||||19940320105800||CARMI|||||19940320105844|||F
- OBX|1|ST|L2565.000^PH, ARTERIAL BLD GAS (R)||7.46||7.35-7.45|H|||F
- OBX|2|ST|L2566.000^PCO2, ARTERIAL BLOOD GAS||28|mm/Hg|35-45|LL|||F
- **NTE**|1|L|BLOOD GAS ANALYSIS CRITICAL VALUE(S) CALLED TO~DR. CARLSON.
- OBX|3|ST|L2567.000^PO2, ARTERIAL BLOOD GAS||83|mm/Hg|80-100|||F
- OBX|4|ST|L2568.000^O2 SAT, ART BLD GAS (R)||96|%|95-99|||F
- OBX|5|ST|L2569.000^BASE EX, ARTERIAL BLD GAS||-2.1|mEq/l|-2.0-2.0|L|||F
- OBX|6|ST|L2570.000^HCO3, ARTERIAL BLD GAS||19.4|mEq/l|22-26|L|||F
- OBX|7|ST|L2571.000^PATIENT TEMP, ABG||96.2|deg F|||F
- OBX|8|ST|L2572.000^MODE, ABG||ROOM AIR|||F

■ <http://www.loinc.org/>

LOINC for coding lab results



- HOME
- BACKGROUND
- RELMA
- DOWNLOADS
- SUBMISSIONS
- NEWS
- MEETINGS
- SLIDESHOWS
- FUNDING SUPPORT
- TERMS OF USE
- LOINC E-MAIL LIST


home

Logical Observation Identifiers Names and Codes (LOINC®)

The purpose of the LOINC database is to facilitate the exchange and pooling of results, such as blood hemoglobin, serum potassium, or vital signs, for clinical care, outcomes management, and research. Currently, most laboratories and other diagnostic services use HL7 to send their results electronically from their reporting systems to

September 29, 2003
Public Clinical LOINC® Tutorial & Committee Meeting
Salt Lake City, UT

Where does LOINC fit in HL7?



Says patient had a serum troponin value of 5 ng/mL, and a serum CK MB value of 24 ug/L.

From www.loinc.org



LOINC's Axes

Table 2. LOINC axes.

Number	Description
1	Component (analyte): e.g., potassium, hemoglobin, hepatitis B antigen
2	Property measured: e.g., a mass concentration, enzyme activity (catalytic rate)
3	Timing: i.e., whether the observation applies to a moment in time or is an average or amount taken over a period of time, as is the case for a 24-h urine sodium concentration
4	System: i.e., type of sample or organ examined: e.g., urine, blood, chest
5	Scale: e.g., whether the measurement is quantitative (a true measurement), ordinal (a ranked set of options), nominal, or narrative (e.g., dictation results from x-rays)
6	Method used to produce the observation, but only when different methods give clinically significant different results




Will 5 axes be adequate?

Table 1. Examples of laboratory LOINC codes and formal LOINC names.

LOINC code	LOINC name (component:property:timing:specimen:scale)
2951-2	SODIUM:SCNC:PT:SER/PLAS:QN
2955-3	SODIUM:SCNC:PT:UR:QN
2956-1	SODIUM:SRAT:24H:UR:QN
2164-2	CREATININE RENAL CLEARANCE:VRAT:24H:UR:QN
1514-9	GLUCOSE^2H POST 100 G GLUCOSE PO:MCNC:PT:SER/PLAS:QN
3665-7	GENTAMICIN^TROUGH:MCNC:PT:SER/PLAS:QN
17863-2	CALCIUM.IONIZED:MCNC:PT:SER/PLAS:QN
2863-9	ALBUMIN:MCNC:PT:SNV:QN:ELECTROPHORESIS

Or will we want to differentiate manual from automated diff results by including the 6th axis, Method?



What will HL7 v3 do to what we do now?

- I doubt that anything will go away
- But I'm not in the HL7 v3 Mood now
- To be continued?



Resources

- HL7 (Health Level Seven) Message Standards for Laboratory Results and Radiology Reports in Ireland 31/7/2003
 - This documents provides ready and succinct description of all the HL7 fields that might appear in a lab results message.
 - <http://www.hebe.ie/ProgrammesProjects/e-Government/SubProjects/GPMessagingProject/Documentation/FiletoUpload,1641,en.pdf>
- LOINC homepage
 - <http://www.loinc.org/>
- McDonald CJ et al. LOINC, a universal standard for identifying laboratory observations: A 5-year update. *Clin Chem* 2003;49:624-33. (access full text via LOINC homepage under "Articles")
- Huff SM et al. Development of the Logical Observation Identifier Names and Codes (LOINC) vocabulary. *JAMIA* 1998;5:276-92. (access full text via PubMed)
- Schadow G et al. Integrating medical information and knowledge in the HL7 RIM. Proc AMIA Symp 2000;764-8.