

caGEDA

UPCI

http://bioinformatics.upmc.edu/GE2/GEDA.html



Why caGEDA?

- Evaluative comparisons of new methods of analysis is rarely conducted and is needed
- Normalization methods are not well understood
- Performance characteristics of tests for identifying differentially expressed genes are understudied
- Optimal combinations of normalization -> feature selection -> sample classification algorithms have not yet been determined
- caGEDA was/is designed with cancer researchers in mind



'Why Not Just Use...'

- BioConductor
- MeV/TM4
- OncoMine
- BRBArray Tools
- GEDP***
- Others...
- Commercial software

- •Please do! Some very nice options!
- •Some require downloads/registration
- •Some require programming
- •Some are not open source
- •Every new microarray data set is another opportunity to identify generally optimized methods of analysis
- *** microarray data repository core to caBIG!!!
- •Training and adoption is efficient with a web application

Normalization	Concept	References	
Reference Gene/ Sample Subset Methods			
Housekeeping Genes	Selection of a set of genes as controls; each value in an array is normalized using the mean of this subset	Lee et al., 2001; Vandesompele et.al., 2002	
'Globalization' Method	Each value in an array is normalized using the global mean of all arrays	Velculescu et al., 1999	
Loess 1: Normalization by self- consistency and local regression	Normalize pairs or groups of arrays relative to each other by iteratively maximizing the consistency of relative expression levels among them. Genes are consistent if their relative expression values do not change after global normalization. The original data are	Kepler et al., 2002	
Iterative Invariant Set Normalization	normalized using the consistent set and local regression Find gene set with unchanged ranks in expression in both groups; use an iterative procedure to identify invariant set as those probes with proportion rank difference (PRD) < 0.003 (low rank) or < 0.007 (high rank genes)	Li & Wong, 2002	
Microarray Sample Pool	Normalize all samples using an ensemble sample (MSP) as the reference array	Yang YH et al., 2002	
Statistical Methods			
Variance Stabilization	Normalization by the arsinh function $h(y) = g \operatorname{arsinh}(a+by)$ with model parameters a and b estimated by likelihood	Huber et al, 2002	
Variance Stabilization	Stabilizes asymptotic variance over the full range of expression intensity. Finds a transformation for a regression model such that the variance is constant over the range of	Durbin et al., 2002	
Dye Channel Control Spot	the dependent variable Expression values normalized by scaling cy5 values so that mean cy5 & cy3 values in control spots are same	Cavalieri et al., 2000	
Loess 2:Local mean normalization	Calculation of local mean (using regression) and distance of this mean from each ratio is the corrected ratio. Results in mean intensity ratio of 1	Colantuoni et al., 2002	
Loess 3:Local variance	Expression ratios made to have same local standard deviation calculated by loess and the intensity is represented as a Z-score	Colantuoni et al., 2002	
Loess 4: Loess Local Regression	Intensity-dependent normalization achieved using the lowess function c(A), specifically $log(R/G)corr = log(R/G)-c(A)$	Yang YH et al, 2002	
Log inverse ratio global normalization	Shift the log ratios by correction factor $log(R/G)$ -correction factor l	Yang YH et al., 2002	
Variance regularization	Normalization factor is calculated using sum of both intensities, which is used to adjust the expression data in its log form	Quackenbush, 2002	
Signal-Dependent Normalization	Center the mean of Cy3 & Cy5 log-ratio distributions	Workman et.al., 2002	
Qspline	Quantiles from target and probe signals used to fit a smoothing B-spline	Workman et.al., 2002	
Spot-Specific Normalization			
Adjustment for slide- specific effect	Ratio-based adjustments: normalize using error factor from simulations; categorical adjustments: use Bartlett's method	Tsodikov et al., 2002	
Spatial Normalization	Subtract local signal estimates from log intensities or log ratios	Workman et.al., 2002	

Test	Reference(s)	
adaptive sign test	Boer et al., 2001	
ANOVA	Kerr et al., 2000; Luo et al., 2002	
BSS/WSS	Dudoit, 2002	
diagnostic metric	Welsh et al., 2001	
discriminative weighting	Bittner et al., 2000	
empirical Bayes method	Newton et al. 2001	
ideal discriminator method	Troyanskaya et al., 2002	
local Bayesian Error test	Baldi and Long, 2001	
log-odds tests	Lonnstedt and Speed, 2002	
neighborhood analysis	Golub et al., 1999	
nonparametric t-test	Garber et al., 2001; Troyanskaya et al., 2002	
perfect discriminator permutation	Park et al. 2001	
Pitman's test	Herwig et al., 2001	
ANOVA with bootstrap variance est.	Black & Doerge, 2002	
significance analysis of microarrays (SAM)	Tusher et al., 2001	
singular value decomposition	Alter et al., 2000;	
	Wall et al. 2001; Ghosh, 2002	
genetic algorithm	Li et al., 2001	
partial least squares	Nguyen and Rocke, 2002	
Welch test	Herwig et al., 2001	
Z-ratio score	Quakenbush, 2002	

Reference(s)
Sultan et al., 2002
Ben-Dor et al., 1999
Quinlan, 1996
Alon et al., 1999
Hastie et al., 2000
Eisen et al., 1998
Kohonen, 1982
Nguyen and Rocke, 2002
Bittner et al., 2000
Shi and Malik
Saitou and Nei, 1987
Li et al., 2001; Theilhaber et al.2002
Bozinov and Rahnenfuhrer, 2002
(e.g., Luo et al., 2002)
Nguyen and Rocke, 2002
Dougherty et al., 2002
Golub et al., 1999; Yeang et al., 2001
Lyons-Weiler et al., 2003





Too many methods.






















Special Capabilities


- Built to facilitate comparisons of methods of analysis via cross-validation + other methods
- Computation validation methods include:
 - Nonparametric bootstrapping
 - Leave-one-out validation
 - Random Resampling Validation
 - − *k*-fold validation (to be added)
 - Efficiency Analysis*** NEW
- Gene Expression Pattern Grid
- Proof-by-Pubmed on the fly

Framework of Evaluation



Credit: Richard Pelikan

Random Resampling




Random data set N1 = N2 = 16; 1100 random 'genes'; t-test


Significance of achieved classification error

Achieved Classification Error, Permutation Achieved Classification, 95th and 99th PACE percentile


Features



Efficiency Analysis




$$O = (2*N3)/(N1+N2)$$




$$O = (2*N3)/(N1+N2)$$

 $N3 = 0$; $O = 0$



$$O = (2*N3)/(N1+N2)$$

 $N3 = N1+N2$; $O = 1.0$



Astrocytoma Progression Markers Early Stage vs. Late Stage

• USA

• Khatua et al.

Early: N = 7

• Late: N = 8

Genes: 8497

Journal: Cancer Res.

CARLES SERVANCES (MILLIAN AND IL 1991)

Overexpression of the EGFR/FKBP12/HIF-2\alpha Pathway Identified in Childhood Astrocytomas by Angiogenesis Gene Profiling^{1,2}

Soumen Khatua, Katia M. Peterson, Kerin M. Brown, Christopher Lawlor, Maria R. Santi, Bunnie LaFleur, Devin Dremman, Dietrick A. Stephan, and Tobey J. MacDonald⁴

Grear for Concer Research Children's Research Restricts [S. E. M. P., C. L. T.) M.). Research Center for Genetic Medicine. Children's Research Durbon (E. M. B., & D., D. a. S.). And Department of Proceedings (M. E. S.). Children's Notices Medicine (Seeding September of Proceedings Medicine Medicine Medicine). The Intellectric Conservity, Medicine, Proceedings (P. M. B.). Seedings (Seedings). Medicine (See

Intense sugingenesis preliferation, a histopathological halbanch distinguidang moligaran from benign astrocytoms, is tital for tumor progrestion. Thus, identifying and carpeling specific patterns; that promote muliganest actrosystems-induced augingments could have raterantial therextracritemes (HGAs) and benign low-grade extracritemes. Hierarchical clustering and petacipal components analysis using only the 150 angiopenetic-related grass distinguished BGA from live-goods astrocytoms in 100% of the complex analyzed, as did nacoparticed analyses using the entire sat of VDR expected gener represented on the arrest, indicating that

in which overall survival comeins less than 30% (2). Thus, novel therapeutic approaches are needed for childhood HQA.

Studies demonstrating the crucial role of angiogenesis in cappet have been a major advance in our understanding of malignant transcr progression (3). One of the key histopathological features that distinspecific beaufit. Expression profiling of 10 childhood surveysums to grade HGA from LGA in intense, increased expression for intense, increase realed that 44 (35%) pases were differentially expressed (37 were ever- with LGA, is associated with increased microvascular density and experied, and 27 were underexpressed) between understand high-grade introductional hypoxia (ii). Then, inhibitors of hypoxia-inducible asgiogenetis factors could be important new thempeutic agents against HOA. In adults, the most commonly described regulators of brain tuzzor-detived suggestatus are VEGF, platelet-derived growth factor. augiopoietin-2, and their respective receptors (5). It is not known to what extent these same regulatory mechanisms exist in pediatric

Germany

van den Boom et al.

Early: N = 8

• Late: N = 8

Genes: 5682

Journal: Am J Pathol.


non-house of Autology, No. 21), No. 3, aptender 2003

Characterization of Gene Expression Profiles Associated with Glioma Progression Using Oligonucleotide-Based Microarray Analysis and Real-Time Reverse Transcription-Polymerase Chain Reaction


Jöre van den Boom," Marietta Wolter," Block Knick * David F. Misak * Andrew S. Youkin, F Daniel S. Wechsler, 1 Clemens Sommer. 9 Guido Reifenberger," and Samir M. Hanashi

Scenarial expression likely plays a cole in astrocytoma. progression. (Am J Partiel 2003, 16):1033-1043)

Diffusely infiltrating astrocytic gliomas are the most common primary beam fumors in adults. 1 These fumors up



Gene Expression Pattern Grid



4b25q06, s1-Hama-sapiens-d0NA-dane-72250-8'-sin Human-mRNA-for-KIAA0026-gene,-complete-cde yc76q06,s1-Homa-sapiens-cDNA-done-21725-3'-sim H. sapiens: UZ1, 1-mRRA vb55h04.s1-Homo-saplens-c0NA-done-75127-3'-sin ym27c08.x1-Homo-sepiena-cDNA-done-49393-3'-sir Human-transforming-growth-factor-beta-induced-ger yi69f12.s1-Homo-sapiens-cDNA-clone-143567-9'-sin Human-gene-for-heterogeneous-nudear-ribonudeo; Homp-sapiens-spliting-factor-(CC1.4)-mRNA,-compl Human-mRNA-for-HMG-1,-complete-cdx yh77hd2.si-Home-sapiens-d0NA-dene-135815-3'-si Human-1-8D-gene-from-interferon-indudble-gene-f Human-spliding-factor-8035-mRNA,-complete-cds ■ Human-mRNA-for-hU1-70K-small-nuclear-RNP-prote Human-26-kDa-cell-surface-protein-TAPA-1-mRNA,-■ Homo-sapiens-mRNA-for-eukaryotic-initiation-factor Human-nudeolar-protein-(822)-mRNA,-complete-cd Human-SPARC/osteonectin-mRNA,-complete-cds yc51a06.s1-Homo-xapienx-cD8A-done-64176-3'-xim Human-mRNA-for-lactoyl-glutathione-lyase yb39b08.x1-Homo-xapiera-d0NA-done-73527-3* H. sapiens-mRNA-for-elongation-factor-2 Human-apolipoprotein-E-(epsilon-4-allele)-gene,-co yd03d04.s1-Homo-sapiens-cDNA-clone-79590-3'-sim ■■■ yp89g11.s1-Hame-sapiens-@NA-dene-194660-3'-si Human-co-beta-glucosidasa-(preactivator)-mRNA,-co Human-bone-morphogenetic-protein-1-(8MP-1)-mR Human-deavage-and-polyadenslation-specificity-fac Human-oANP-response-element-regulatory-protein-Human-mRNA-for-lactate-dehydrogenase-A-(LDH-A, gapdh-5 4h87e06, s1-Homo-sepiens-d0NA-done-136738-3'-si yf84g10.s1-Hemo-sapiens-d00A-done-29328-2 yq67e02.s1-Home-sapiens-d0NA-dene-200858-31 Human-lyeazyme-mRNA,-complete-cde Human-myelin-associated-glycoprotein-(MAG)-mRNA Human-deavage-and-polyadenylation-specificity-fac yj03f03.s1-Homo-sapiens-cDNA-clone-147677-8'-sin yi48e06.x1-Homo-xepienx-cDNA-clone-142498-3'-xi ya88f84, s1-Hemo-sapiens-dDNA-done-68767-3 Human-mucin-2-(MUC2)-mRNA-sequence gu67d10, g1-Home-sapiens-d0NA-dane-228967-2'-si Human-Intestinal-trefoil-factor-(HITF)-gene,-exon-8 yj03f03.x1:Homo:xepienx:cD8A-clone-147677-3'-sim yb24di1.si-Home-sapiers-d0NA-dene-72117-8' yc22h05.s1-Homotxapiens-cDNA-done-81465-3'-xin b-actin-3 yb27c07,s1-Homorsepiens-cDNA-clone-72396-3'-s Human-CO-029 yu75h09.s1-Homo-sapiens-d0NA-done-239681-3'-si Home-capiene-colon-mucosa-acsociated-(DRA)-mRI ub87d09.s1-Homo-sapiens-d0NA-done-78161-8'-sin uf18=03.x1-Homo-sapiens-cDRA-done-127228-3'-sir vh88b10.s1-Home-saplens-cDNA-done-136795-3'-: yb25f04.x1-Homoraspiena-dDNA-done-72223-3'-xim Human-ribosomal-protein-825-mRNA,-complete-ods yb22e07, s1-Home-sepiens-d0NA-dene-71940-9'-sin ye38e05, s1-Home-sapiens-d0NA-dane-120032-2'-si yu67a12.s1-Homo-sapiens-dDNA-done-238846-31 ya91a12, z2-Home-sapiens-d)NA-dene-69022-3'-sin H. sapiens-p27-mRNA ub+0c08,x1-Homo-xapienx-cD8A-done-73646-3'-xim Home-sapiens-guanylin-mRNA,-complete-cds Human-carbonic-anhydraxe-IV-mR8A,-complete-cdx yg31b10.s1-Home-sapiens-d0NA-dene-197371-3'-si yc48a10.s1-Homo-sepiens-cDNA-done-83898-3'-sim yh56 a10, x1-Home-capiens-d0NA-dene-122722-2'-si ye49f12.s1-Homo-sapiens-dbNA-done-121108-3'-sir Human-cysteine-rich-protein-(CRP)-gene, -exons-5-Human-cysteine-rich-protein-(CRP)-gene,-exons-5-Human-smooth-musda-cell-celponin-mENA,-compleyd72b02.si-Homo-saplens-cDNA-done-113739-3'-si Human-cysteine-rich-protein-(CRP)-gene,-exons-5-; yo70c01.s1-Homo-sapiens-cDNA-done-183264-3'-st yb43f08.s1-Homo-sepiens-dDNA-done-73959-S'-sim yc50dQ4.s1-Homo-sapiens-cDNA-done-94192-2'-cim yn63h10.s1-Home-sapiens-d0NA-done-178155-3'-si yb57d09.x1-Home-xapienx-d0NA-done-75161-3'-xin H. sapiens-ACTB-mRNA-for-mutant-beta-actin-(beta'-Over in A. Under in B Over in A Under in A Under in B Over in B Under in A. Over in B. G. Others

Not Expressed


'G' = other

Some tests lead to more sizeable 'G' group which, while statistically significant, exhibit no coherent signs of differential expression in most samples.

Outliers or conflicting patterns of differential expression.

(colon cancer data set, t-test, cut-point = 4.0)







- Enhance!
- Integrate and Interoperate!
- Annotate!
- Blow it up!
- Characterize and represent
 - Data models
 - Schema
 - UML Diagrams:
 - Use case diagrams
 - Activity diagrams
 - Sequence diagrams
 - Package diagrams...



Enhance!

- Increase data format diversity tolerance
- Add outlier spot detection, adopt existing QC criteria
- Add normalization (e.g., DWD), tests, classification methods
- Apply Jprogram (Duke) to allow assimilation of R projects
- Add pathway analysis and interaction analysis capabilities (cMAP, cPATH, cytoScape...)



Integrate and Interoperate!


- SPOT, SPROC, LIMS projects, OncoMine, FDGP could produce data dumps in caGEDA formats or adopt html interface that finds an active caGEDA server (local or on the grid) for on-the-fly analysis
- caGEDA could output in formats or make direct calls to:
 - GoMiner
 - cPATH
 - GKB (Reactome project)



Annotate!

- Five components:
 - English text description
 - Mathematical description
 - Pseudocode
 - Source code
 - Related literature

Blowing it up...





Demo

• http://bioinformatics.upmc.edu/GE2/GEDA.html


Development (credit: S. Patel)

caGEDA application development is an iterative software development approach that leverages elements from Rational Unified Process (RUP). Use cases for the application are developed using the expertise available at the research center. Once the use case analysis is completed, an iterative functional design and development process is applied, which allows for rapid and segmented development of the application. During the iteration, all the software development activities are executed. The artifacts associated with each functional iteration includes: detailed use cases describing the function; class and sequence diagrams; a system architecture diagram; the actual software code; a project plan describing subsequent iterations; and a test plan for software validation.

• UML modeling and use case development is performed using UML modeling tool from Rational Rose. Source code is developed using the Java programming libraries for Servlet, JSP and EJB. We use Apache software's Ant to assist the software build process. All the server side software components are tested on JBoss application server. All the software components used in development the GEDA application are freely available on the Internet.

Application Architecture

caGEDA conforms to n-tier architectural design that include several layers. A presentation layer includes a web application server that transforms the request coming from the Internet browser in to the calls to the business logic and provides programmatic access to the application. A Business layer can communicate to the standalone application client directly using RMI-IIOP protocol or using CORBA. Presentation layer objects communicate with the business layer objects using RMI-IIOP protocol. Since the communication with the business layer can be done using CORBA even a non-java application client can make use of the services provided by caGEDA.



Presentation tier involves one or more web servers, each responsible for interacting with end user. The presentation tier displays the requested information in HTML to the end user; it also reads and interprets the user's selection and makes invocations to the business tier's components. The implementation of presentation tier uses Servlets and JSPs

Business tier consists of multiple EJB components running under the hood of EJB container/server. These are reusable components that are independent of any user interface logic. We should be able to, for example, take our business tier and port it to different presentation tier (such as application client) with no modification. Our business tier is made up of session, entity and messagedriven beans.

Data tier is where the permanent data resides. With use of entity beans, we can leave our options to use virtually any database of choice. Switching to database of a particular choice should be seamlessly achievable.

Software Life Cycle: Iterative approach

High level requirement analysis

Scope

Data Format

Data Pre-processing

Feature Selection

Prediction

Computational Validation

Data Visualization

Databases

Security

Architecture and Design

Estimation & Schedule

Iterative Design & Development

Testing

Continuous testing and integration

Coding standards

Implementation

Testing

Beta Testing


Feedback

Deployment

Deployment plan

User documentation

Bug reporting/tracking system



Deployment